










# SIDE RAIL QUICK START GUIDE


## INCLUDED



**WARNING**  
DO NOT TILT or LEAN base on legs when turning base over.



# SIDE RAIL QUICK START GUIDE



**Zone 1: Under the Rail** - The gap under the rail between the mattress, may allow for dangerous head entrapment. The FDA recommended space is less than 4.75 in (12 cm).

**Zone 2: Between the Rail and the Mattress** - This area is the space between the inside surface of the bed rail and the mattress, and if too big it can cause a risk of head entrapment. The FDA recommended space is less than 4,75 in (12 cm).

**Zone 3: Between the End of the Rail and the Side Edge of the Head or Foot Board** - A gap of less than 12.5 in (32 cm) between the end of the bed rail and the side edge of the headboard or footboard can present the risk of entrapment of the head, neck or chest.

## Warnings

There is a risk of entrapment associated with all bed rails and similar bedside mobility aids. Entrapment occurs when an individual becomes wedged between a bed rail/handle and another object, such as the mattress, potentially resulting in serious injury or death.

## Before Installation:

- Read the enclosed **Entrapment Prevention Guidelines** before installing or using this product. For further information, visit <https://www.cpsc.gov/safety-education/safety-guides/furniture-furnishings-and-decorations/adult-bed-rails> and [www.fda.gov/bedsafety](http://www.fda.gov/bedsafety)
- Ensure that it is safe to use on your bed and mattress. This product can only be used on beds and mattresses that comply with the specifications in the enclosed **Mattress Conditions & Safety Components** and **Assembly Instructions**.
- Do not use this product with beds that may cause a gap between the mattress and the bed rail, such as waterbeds, lightweight foam, air, soft pillow-top mattresses, or other very soft mattresses.
- Review all specifications and instructions prior to beginning installation.
- Install this product per the enclosed **Entrapment Prevention Guidelines**. If you have any questions about whether the product is safe for use on your bed or how to install or use it, please call us at 314-499-8869.

## Usage Purpose:

- This product is intended to assist disabled persons by providing minimal to moderate weight support while standing from bed.
- This product is not intended to be used as a device to restrain individuals from getting into and out of bed.



# SIDE RAIL QUICK START GUIDE

## Additional Warnings:

(a) This product is equipped with components to securely attach it to a Symphony Sleep adjustable frame ONLY. It is not compatible with any other bed frame or adjustable bed. Failure to properly install and use these components significantly increases the risk of entrapment and may result in personal injury or death. NEVER remove these components or use the product without these components installed as per the enclosed **Assembly Instructions**.

(b) This product must be installed on the adjustable bed. Never use this product without a Symphony Sleep adjustable bed.

(c) Ensure there is no gap between the product and the side of the mattress before use. A gap in this area means improper installation. Stop using immediately if the product shifts out of its original position, a gap is noticed between the side of the mattress and the product, or if the product is damaged, broken, or missing parts.

(d) This product is NOT intended for use in a nursing home, assisted living center, or similar facility.

(e) This product is for adult use only. It is not intended for use by individuals who have paralysis, symptoms of dementia, sleeping disorders, incontinence, pain, uncontrolled body movement, or those unable to get out of bed and walk safely without assistance. It should not be used by individuals who suffer from confusion, restlessness, terminal restlessness, or those under the influence of medications, drugs, or any substance that could impair balance or judgment, or any other unforeseeable reasons that could affect the user's physical and mental ability to safely use this product. If you have questions, consult a physician before using this product.

(f) When changing bed linens (bed covers, blankets, sheets, etc.), always confirm that the product is installed correctly and re-tighten the clamps to eliminate a gap between the product and the side of the mattress.

(g) If the product is removed and reinstalled, follow the enclosed **Assembly Instructions**.

## Individuals at Higher Risk:

- People with Alzheimer's disease or dementia, or those who are sedated, confused, or frail are at increased risk of entrapment and strangulation.
- Individuals attempting to climb over the product are at increased risk of injury or death from falls.

## Do Not Use:

- NEVER use unless the product is tight against the mattress without gaps and at least 12 ½" (32 cm) from the headboard and footboard.
- NEVER use with children.
- NEVER use for toddlers or infants.